

Supply and installation of specially designed Laboratory Furniture.

The manufacture and installation of the furniture involved the plinths of the furniture being laid as a separate item, to enable the flooring to be covered up the plinth, thus ensuring a completely sealed floor.

Following installation of the flooring, the base units were then attached to the pre-laid plinths and VELSTONE work surfaces installed. The use of VELSTONE enabled the client to have completely joint free work surfaces, which in this instance was of paramount importance.



Industrial Research Facility Bespoke Furniture

InterFocus



lab line